

# NEW KU-3™ Viscometer

for Paints, Coatings, and Inks

**ASTM D562 Compatible**  
(industry specification)

**New Magnetic Spindle Coupling**

**New User Interface**  
with touch control & auto-start

**New Single Piece Can Adapter**  
for pint and 1/2 pint cans  
*Accommodates quart cans*

**Easy to use**  
no weights, simplifies an  
established test procedure

**LED Display Info:**

- Krebs Units
- Gram Units  
(Weight)
- Centipoise\*

Select Krebs or Grams or Centipoise

**Lock-In Test Results**  
with Hold button

**Accuracy: ±1.0% of range**

**Repeatability: ±0.5%**

**Standard Krebs Spindle**

**Measurement range:**  
40 to 141 KU, 32 to 1099 gm,  
and 27 to 5274 cP\*



## What's Included?

- Instrument
- Krebs-type Spindle (p48)
- Adapter for Quart, Pint and Half-Pint Cans

## Optional Accessories

- Paste Spindle (p48)

## Applications

- Paints
- Coatings
- Adhesives
- Inks
- Pastes

\*Centipoise values based on the conversion from Krebs Units as defined in the ASTM D562.